



Human genomic control DNA

Mutation in α -globin gene (α -thalassemia)

CERTIFICATE OF ANALYSIS

Product Information

Product name :	Normal Control DNA for α -Thalassemia
Cat. No. :	K5411
Concentration:	300 ng/ μ l
Storage condition :	-20°C
Mutation type:	$\alpha\alpha$ IVSI-4(A>G) / $\alpha\alpha$

QC Results

Test Items	Specifications	Results
Absorbance ratio	A260/A280>1.8 A260/A230>2.0	Passed
Electrophoresis	Average size larger than 50 kb	Passed
Sequencing	No other mutation in the gene	passed
Enzyme activity	<i>Eco</i> RI digestion	passed

Caution:

This product is pure human genomic DNA, and must be kept away from PCR products to avoid contamination

Sequencing result :

